



THE DENTAL FACTORY

Passive, Active and Functional
Fixed/Removable Orthodontic appliances.

Splints / Orthotics

Teeth Alignment trays


(We will work with you through multi-stage cases)

Sports Mouthguards

Catering for all sports

Any custom design you may have... We love a challenge!

Time Saver PRE-Booking.. call us on **1300 44 15 84** to pre-book any future cases, and in doing so, save turnaround time by a considerable amount!





BIOBLOC STAGE I

Adjustments


PLIERS

- Adams
- Three point








RETENTION

- Adams pliers
- Sub-gingival
- Keep length for future adjustments







DURING EXPANSION



ADJUSTING POINTS

- Utilise the positioning of the pre-bent curves in the wires in order to activate them
- Always make sure you have enough retention before activation
- Maintain a guide plain using the labial Archwire as a possible reference





ADJUSTING POINTS